





Regenerating Vision with Light

Infectious keratitis and **corneal ulcer** are common diseases in dogs, cats, horses and birds.

They can be rapidly complicated by *corneal melting* leading to visual loss.



dog



horse



bird

cat

Corneal photochemotherapy with **VETUVIR**

is the **new antimicrobial therapy** for efficiently treating **infectious ulcerative keratitis** with **riboflavin** and **UV-A light**.



Indication for use: Infectious keratitis and septic corneal ulcer in dogs, cats, horses and birds.

How does it work?

Corneal photochemoterapy is performed by irradiating the cornea, which has been soaked with **VISIOFLAVIN** riboflavin solution for 20 minutes, using **VETUVIR**, the advanced UV-A device.


- 1) Soak the cornea with VISIOFLAVIN
- 2) Irradiate the cornea with VETUVIR

VISIOFLAVIN



VETUVIR

VETUVIR and **VISIOFLAVIN** are for veterinary use only. Trademarks registered by Vision Engineering Italy srl.



Unique features

Two UV-A irradiance settings: 10 mW/cm² and 30 mW/cm²

Fully portable

Easy-to-use

Bluetooth remote control

*Efficiently used either
in the operating room
or in the stable*



Manufactured according to Dir. CEE 2001/95, 2006/42/CEE and 2014/30/UE

Riboflavin with enhancers

0.1% riboflavin solution

- *hypotonic*
- *dextran-free*
- *with enhancers*



Clinical cases

Courtesy Prof. R. Gialletti

Horse with chronic corneal abscess
(pseudomonas+ and streptococcus+)
Resistant to antibiotic therapy.



1 week after treatment:
Almost complete resolution
of the abscess.





5 weeks after treatment:
(pseudomonas- and streptococcus-)
Vision is restored.

Clinical cases

Courtesy Dr. C. Giordano

Corneal photochemotherapy

Persian cat
with ulcerative keratitis



1 month after treatment



2 months after treatment
The left cornea is almost CLEAR

Clinical cases

Courtesy Prof. C. Peruccio

French bulldog with corneal melting and descemetocoele.
Resistant to antibiotic therapy.



2 weeks after treatment:
The cornea becomes
more transparent.



Photochemoterapy



1 day after treatment:
(no more descemetocoele)



4 weeks after treatment:
Visual function is preserved.



ORDERING AND INFORMATION

Vision Engineering Italy srl

Phone: +39 331 34 38300

Email: visioengvet@visioeng.com

www.vetuvir.com

VISIOFLAVIN



IMPROVING EYE CARE IN VETERINARY MEDICINE

